## ***MANIPOLO DA LABORATORIO***

CERTIFICATO DI GARANZIA  
GUARANTEE CERTIFICATE  
CERTIFICAT DE GARANTIE  
GARANTIEZERTIFIKAT  
CERTIFICADO DE GARANTIA

L'ODONTOTECNOLOGIA

FARO s.p.a  
Via Faro 15  
20876 Ornago (MB) Italia  
Tel: 039/6878.1 (r.a.)  
Fax: 039/6010540

Internet: <http://www.faro.it>

## ***MANIPOLO DA LABORATORIO***

MANUALE DI ISTRUZIONI  
INSTRUCTION MANUAL  
MODE D'EMPLOI  
GEBRAUCHS ANWEISUNGEN  
MANUAL DE MANEJO

L'ODONTOTECNOLOGIA

**ED. 9/11**  
**Rev. 0**  
cod. 988.060

## 6 mesi - months - mois - monaten - meses

*nome - name - nom - vorname - nombre*                      *cognome - surname - prenom - nachname - apellido*

*indirizzo - address - adresse - auschrift - direccion*                      *città - town - ville - ort - ciudad*

*data di produzione - production date*                      *data d'acquisto - purchase date - date d'achat - einkau-  
fdatum - fecha de compra*

Manipolo/Handpiece:

- 010.000       012.000
- 011.000       014.000
- 011.020

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---

12

**FARO laboratory handpieces conform to  
CEE machine directive 89 / 392**

## **IMPORTANT !**

**Before installing the handpiece read this manual carefully**

13

## **TABLE OF CONTENTS**

Safety rules	page	15
Features	page	17
Operating Instructions	page	18
Maintenance	page	21
Troubleshooting	page	22
Technical particulars	page	23
Guarantee Certificate	page	24

14

## **SAFETY RULES**

- Follow all the instructions set forth in this manual.
- ◆ To avoid undue wear on the handpiece do not brake rotation by using the collet opening lever for any reason.
  - ◆ Do not open the collet when the shaft is rotating.
  - ◆ Do not use the handpiece with the front ring nut partially unscrewed.
  - ◆ The handpiece must not rotate without the fitted tool.
  - ◆ It is recommended to wear protective clothing, eyeglasses and gloves while using thehandpiece.
  - ◆ Do not attempt any maintenance while the shaft is rotating.
  - ◆ Insert only tools of the correct diameter in the handpiece collet.

15

- ◆ Do not attempt to repair the handpiece personally but only perform those operations explained in the manual; in all other cases call on our Service Department.
- ◆ Use only tools designed to resist a rotation speed higher than that of the handpiece (see the tool manual).
- ◆ When working on crumbly materials it is recommended to use an exhaust located near the place the tool is working in.

16

---

## FEATURES

- ◆ Rotation axle on self-adjusting bearings.
- ◆ Easy tool change.
- ◆ The handpiece is available in the following models:

handpiece slip connection with collet $\varnothing$ 2,35 mm	(cod. 010.000)
handpiece flexible connection with collet $\varnothing$ 2,35 mm	(cod. 011.000)
handpiece duplex connection with collet $\varnothing$ 2,35 mm	(cod. 011.020)
handpiece USA connection with collet $\varnothing$ 2,35 mm	(cod. 012.000)
handpiece slip connection with collet $\varnothing$ 3 mm	(cod. 014.000)

---

17

---

## OPERATING INSTRUCTIONS

Check packing to see that all parts of the laboratory handpiece are received, as follows:

- ◆ Handpiece
- ◆ Annexed documentation

### INSTALLATION OF HANDPIECE

- ◆ **Fit the handpiece by pushing it into the flexible shaft connection by taking care to position the reference notch with the connection system.**


---

18

---

### CHANGING TOOL

To insert the tool shank in the handpiece collet move the lever (2) to position B (collet open). To lock the tool move the lever (2) to position A (collet close).

**If you have difficulty in withdrawing the tool shank, do not force it but on the contrary push it inward to facilitate opening of the collet.**

---

19

---

**Run the handpiece without load for a few seconds to make sure it rotates properly.**

**HANDPIECE MUST ROTATE ONLY WITH THE FITTED TOOL.  
NEVER BRAKE ROTATION BY USING THE LEVER.**

---

20

---

## MAINTENANCE

Periodically loosen the front knurled ringnut to bare the collet locking system.  
Move the collet locking lever to position B (fig. 1) and blow a strong air jet into the collet.

### SELF-LUBRICATING HANDPIECE

**Do not lubricate**

---

21

---

## TROUBLESHOOTING

**The tool shank does not lock in the closed collet.**

- Make sure the shank diameter is compatible with the collet.
- Call our Service Department to have the collet calibrated.

**The tool does not rotate correctly and the handpiece tends to heat up.**

- Clean the handpiece as explained under MAINTENANCE.
- Call our Service Department.

---

22

---

## TECHNICAL PARTICULARS

Maximum rotation speed	: 22000 rpm
Handpieces various models	: for tools diameter 2,35 / 3 mm
Maximum transmissible torque	: 4,9 Ncm
Noise level	: < 60 dB (A)
Vibrations	: < 2,5 m / s <sup>2</sup>
Connection transmission shaft side,	: type slip / flexible / duplex / USA
Dimensions	: length 154 mm - diameter 17 mm (mod. 010.000)
Weight	: 110 grams (mod. 010.000)

---

23

---

## **GUARANTEE CERTIFICATE**

FARO offers the final customer a **6 month** guarantee starting from the date of purchase.  
Repairs under guarantee must be performed at FARO; expenses and transport risks are at the risk of the purchaser.  
**Repair under guarantee is considered valid only when:**

- **all sections of the certificate have been filled in and sent in advance to FARO by Fax (039.6010540).**

The guarantee covers faults due to the bad quality of the material or manufacturing defects; in the case of valid claims, the guarantee covers free repair or replacement. **Claims for damages and/or interest are excluded.**  
The guarantee is not considered valid, at the sole discretion of FARO, if the fault is due to tampering, damage, incorrect use, improper maintenance and normal wear and tear.

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