D GARANTIEZERTIFIKAT

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GB GUARANTEE CERTIFICATE

Faro offers the final customer a **6 month** guarantee starting from the date of purchase.

Repairs under guarantee must be performed at FARO; expenses and transport risks are at the risk of the purchaser.

Repair under guarantee is considered valid only when:

- **all sections of the certificate have been filled in and sent in advance to FARO by Fax (039.6010540).**

The guarantee covers faults due to the bad quality of the material or manufacturing defects; in the case of valid claims, the guarantee covers free repair or replacement. **Claims for damages and/or interest are excluded.** The guarantee is not considered valid, at the sole discretion of FARO, if the fault is due to tampering, damage, incorrect use, improper maintenance and normal wear and tear.

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**The laboratory motor 22.000 r.p.m.
complies with the directives
89/336 CEE EMC 73/23 CEE DBT**

IMPORTANT!

Before installing the motor, read this manual carefully.

Modello

Laboratory motor 22.000 r.p.m. with electronic control

TABLE OF CONTENTS

Safety Rules	page	10
Description of Parts and Controls	page	11
Operating Instructions	page	12
Maintenance	page	13
Troubleshooting	page	14
Technical Particulars	page	15

SAFETY RULES

Laboratory motor with electronic control and flexible shaft conceives for the following use:

Dental laboratories, jewellery laboratories and mechanical tools

- Follow all the rules set forth in this manual.**
- The use of this motor must be carried out by well trained personnel only.
- During maintenance operations disconnect the power cord from the power outlet.
- Do not perform maintenance on the flexible shaft with the motor turning.
- Do not allow foreign objects in the flexible shaft connector.
- When using the motor it is recommended to wear appropriate clothing, eyeglasses and protective gloves.
- Screw down tightly the ring nut fastening the flexible shaft to the motor.
- During use avoid blocking the burr.
- Use only handpieces capable of withstanding the top rotation speed of the motor.
- Use only handpieces with connection adapted to the flexible shaft.
- Do not perform maintenance operation different from those set forth in this manual. Any work not contemplated by the documentation could involve risk. For this type of operation call on the Service Department.

**Never attempt maintenance on motor
with power on:**

DANGER OF ELECTROCUTION!

DESCRIPTION OF PARTS AND CONTROLS

1. Motor
2. Motor Power Cord
3. Flexible Shaft
4. Rheostat
5. Mains Power Cord

fig. 1

OPERATING INSTRUCTIONS

Check the packing to make sure all the components of the Suspended Motor are received. They are:

- Rheostat
- Motor
- Flexible Shaft
- Annexed Documentation

INSTALLATION OF MOTOR

- Secure the clamp stand (optional) to the bench firmly and safely.
- Secure the motor to the clamp stand locking it with the screw pin provided.
- Connect the Flexible Shaft to the Motor and the Handpiece to the Flexible Shaft.
- Plug the cord into the mains power outlet.

Before plugging in the power cord make sure the mains voltage is the same as that shown on the rating plate of the motor.

ADJUSTING ROTATION SPEED

- Select speed by moving the rheostat lever to the right for slow and to the left for fast (fig. 1, pag. 11).

STARTING THE CONTROL UNIT

- Move the rheostat lever from right to left (fig. 1, pag. 11).

Run motor without load for a few seconds to make sure it operates correctly.


MAINTENANCE

REPLACING BRUSHES

When it becomes necessary to replace the brushes proceed as follows:

- Disconnect power cord from outlet.
- Remove top cover using an Allen wrench to extract the two screws holding it (fig. 2).
- Withdraw the tongue "A" and replace the brushes (fig. 2).
- Assemble the top cover with the two screws (fig. 2).

**BE CAREFUL NOT TO REVERSE THE BRUSH WIRES.
THIS WOULD CAUSE THE MOTOR TO ROTATE IN
THE OPPOSITE DIRECTION.**

REPLACING FUSE

- The fuse is inside the rheostat. To reach it remove the hexagonal screw on the top of the cover.

TROUBLESHOOTING

Motor does not turn.

- Make sure power cord is plugged into outlet and that there is power at the outlet.
- Check condition of fuse located in rheostat.
- Disconnect flexible shaft from motor and try again.
If the motor runs now, the trouble is in the flexible shaft or the handpiece.
- Check brushes for wear and change if necessary (see page 13).
- Call our Service Department.

Motor rotates correctly but power is not transmitted to flexible shaft.

- Disconnect flexible shaft from motor and check condition of transmission gears between shaft and motor.
- If replacement appears necessary call our Service Department.

Power is transmitted normally to flexible shaft but the handpiece axle does not turn.

- Disconnect handpiece from flexible shaft and check condition of transmission.
- If replacement appears necessary call our Service Department.

TECHNICAL PARTICULARS

SUSPENDED MOTOR

Rotation Speed	: 22.000 rpm
Voltage input	: 230 V ac (or 115-120-240 V ac)
Mains frequency	: 50/60 Hz
Power	: 300 W
Maximum absorbed current	: 3 A (for 230-240 V) 6 A (for 115-120 V)
Maximum torque	: 12 N·cm
Noise level	: < 60 dB (A)
Vibrations	: < 2,5 m/s ²
Weight of motor	: 4 kg (8.8 lb.)
Class	: II (double insulation)
Protection grade	: common
Operation type	: continuous with alternating load
Working temperature	: 0÷40°C (32,104°F)

RHEOSTAT

Rotation speed	: 8.000÷22.000 rpm
Power supply	: See motor specifications
Fuse for 230 or 240V	: T 3.15A
for 115 or 120V	: T 6.3A

RHEOSTAT "CONTRONIC"

Rotation speed	: 3.000÷20.000 rpm
Power supply	: See motor specifications
Fuse for 230 or 240V	: T 3.15A

FLEXIBLE SHAFT

Rotation speed	: 22.000 rpm
Power supply	: See motor specifications
Coupling type	: towards motor ⇒ ring nut towards handpiece ⇒ normal slip

Disposal container. Dispose the device observing the standards about the differentiated collection of electrical material.

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