01

# MP 3020

MANUALE D'USO  
INSTRUCTION MANUAL  
MODE D'EMPLOI  
GEBRAUCHSANLEITUNG  
MANUAL DE USO

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## TABLE OF CONTENT

|     |   |    |
|-----|---|----|
| 1   | GENERAL SAFETY RULES .....                            | 14 |
| 1.1 | Intended Use.....                                     | 14 |
| 1.2 | Installation.....                                     | 14 |
| 1.3 | Maintenance and repair .....                          | 14 |
| 1.4 | Contra Indications .....                              | 15 |
| 1.5 | Cleaning and disinfections .....                      | 15 |
| 1.6 | Disposal of the device .....                          | 15 |
| 1.7 | Storage Conditions .....                              | 15 |
| 2   | DESCRIPTION OF THE DEVICE .....                       | 18 |
| 2.1 | Description of functioning.....                       | 18 |
| 2.2 | Description of the parts .....                        | 18 |
| 3   | INSTALLATION AND CONNECTIONS.....                     | 19 |
| 3.1 | Assembly of the lamp on the dental unit.....          | 19 |
| 3.2 | Power connections.....                                | 19 |
| 3.3 | Remote cable connections .....                        | 19 |
| 4   | INSTRUCTIONS FOR USE .....                            | 20 |
| 4.1 | Switching on/off adjustment.....                      | 20 |
| 4.2 | Adjustment .....                                      | 20 |
| 4.3 | Lamp With Remote Control .....                        | 20 |
| 4.4 | Overall movements.....                                | 21 |
| 4.5 | Acoustic Warning Signals .....                        | 21 |
| 5   | MAINTENANCE AND PERIODIC CHECKS.....                  | 21 |
| 6   | CLEANING, DISINFECTION AND STERILIZATION.....         | 22 |
| 6.1 | Chemical products not allowed .....                   | 22 |
| 6.2 | Cleaning the mirrors.....                             | 22 |
| 6.3 | Cleaning of handles, arms and joints.....             | 22 |
| 6.4 | Sterilization Of The Handles .....                    | 22 |
| 7   | TROUBLESHOOTING.....                                  | 23 |
| 8   | TECHNICAL SPECIFICATION.....                          | 23 |
| 8.1 | Optical performances at 1000 mm from the target ..... | 23 |

## SYMBOLS

Dispose of the device in compliance with provisions for separate collection of electrical material.

Device is marked Class II of insulation.

Consult Instruction for use.

**CAUTION**

## **1. GENERAL SAFETY RULES**

### **1.1 Intended Use**

MP3020 lamp has the function of lighting anatomic parts of the human body in following fields:

- Podology,
- Cosmetics,
- Veterinary,

The lamp is not designed to be used during surgical intervention.

The articulated arm and hinge joints of the head make it possible to position the light beam according to user's needs.

The lamp must only be used by qualified personnel after the reading of this manual.

### **1.2 Installation**

Verify that the power supply voltage, shown on the label, is the same as the mains voltage. The device must be installed only by qualified personnel, competent in installation of electrical equipment in medical rooms.

The dental light must be installed on a specific control and supply device, such as a dental unit, or with an electrical system that complies with standard IEC 364-1 and "national regulations for installation of electrical systems in areas for medical use".

The device must be installed with an omnipolar shut-off switch that disconnects it from the mains and is conformant to IEC/EN 61058. It is necessary to install a green signal light to indicate when the lamp is switched on.

The device is not suitable for installation in environments containing flammable gas or rich in oxygen.

### **1.3 Maintenance and repair**

Do not perform any maintenance operation on the lamp when it is connected to the power supply or while the patient is present.

Disconnect the power chord from the mains before starting maintenance procedures.

The user is not allowed to perform maintenance operations not described in this manual. Any operation not indicated in the manual can compromise the safety characteristics of the device.

For the maintenance or repair of this device use only Faro's original spare parts.

Do not insert any objects or tools which could touch parts under voltage into the slots on the lamp head.

### **1.4 Contra Indications**

Do not address the light beam directly in the eye.

Following patients can be particularly subject to photobiological risk: children and adults with eye diseases, people using photosensitive substances (drugs, cosmetics, etc) that can release phototoxic substances.

These patients must not stare at the light beam and must use appropriate protection devices and precautions.

Do not leave small components of the device in unattended or within reach of people at risk (children) because they could be a source of danger.

The high electromagnetic energy of the electro surgery unit can affect the functioning of the other Medical Equipment. Make reference to the instruction of the Manufacturer of such equipment for the recommended distances that must be held from other devices.

## 1.5 Cleaning and disinfections

For cleaning and disinfection DO NOT USE detergents-disinfectants containing:

- Ammonium Hydroxide;
- Sodium Hydroxide;
- Methylene chloride;
- Methilalcohol;
- Abrasive substances;
- Trichloroethylene;
- Benzene;
- Turpentine;

Failure to comply with this indication might cause:

- risk of breakage or damage of plastic parts
- risk of damage of the mirrors

See dedicated section for the allowed products.

In case of doubt, please contact FARO's Service.

Do not spray any liquid directly on the back of the head of the lamp. It could lead to a shortcircuit

of electric-electronic parts. To prevent cross contamination between patients follow instructions in section 6 of this manual.

## 1.6 Disposal of the device

The lamp is a Electrical Electronic device.

The scrap or disposal of this device must comply with national applicable regulations.

If necessary, contact specialised companies that have the necessary authorisation for the transportation and the scrapping of Electrical Electronic devices.

Separate the materials based on their typology (ex. metals, rubber, plastic, etc.).

## 1.7 Storage Conditions

The lamp's packaging is suitable to protect it correctly from penetration of external agents.

The device in its original packaging can be transported or kept in storage for a period of 15 weeks under the following environmental conditions:

- Room temperature between 20°C and 70°C
- Relative humidity between 10% and 90%
- Atmospheric pressure between 500 and 1060 mBar.

The device must be used under the following atmospheric conditions:

- Temperature between 10°C and 40°C
- Relative humidity between 30 and 75%
- Atmospheric pressure between 700 and 1060mbar

## REQUIREMENTS FOR ELECTROMAGNETIC COMPATIBILITY

This medical device requires particular precautionary measures to ensure electromagnetic compatibility, and shall be installed and used in compliance with information provided in the accompanying documentation.

| Manufacturer's guidelines and statement – Electromagnetic emissions  |            |   |
|--|------------|---|
| The lamp MP 3020 is designed to function in the electromagnetic environment specified below. The client or user must ensure its use in the said environment. |            |   |
| Emission tests   | Compliance | Electromagnetic environment - Guidelines  |
| RF Emission CISPR15  | Compliant  | The lamp MP 3020 uses RF energy only for its internal function. Therefore its RF emissions are very low and most likely do not cause any interference in neighbouring electronic devices. |
| RF Emission CISPR15  | Compliant  | The lamp MP 3020 is fit for use in all buildings, including domestic ones and those directly connected to the public low voltage supply network that feeds buildings for domestic use.    |
| Harmonic emission  | Class C    |   |
| Voltage fluctuations/flicker emission  | Compliant  |   |

## ELECTROMAGNETIC IMMUNITY

| Manufacturer's guidelines and statement – Electromagnetic immunity   |  |  |
|--|--|--|
| The lamp MP 3020 is designed to function in the electromagnetic environment specified below. The client or user must ensure its use in the said environment. |  |  |
| Immunity test  | Compliance   | Electromagnetic environment - Guidelines   |
| Electrostatic discharge (ESD)<br>IEC/EN61000-4-2   | ± 6kV contact<br>± 8kV air   | The floor must be in wood, concrete or ceramic. If the floor is covered with synthetic material, relative humidity should be at least 30%.   |
| Electrical fast transient/burst<br>IEC/EN61000-4-4   | ± 2kV power supply<br>± 1kV for input/output lines                                       | The quality of supply network voltage should be typical of commercial or hospital environments.  |
| Surge<br>IEC/EN61000-4-5   | ± 1kV differential mode<br>± 2kV common mode   | The quality of supply network voltage should be typical of commercial or hospital environments.  |
| Voltage dips, short interruption and voltage variation<br>IEC/EN61000-4-11   | < 5% Ut for 0,5 cycle<br>40% Ut for 05 cycle<br>70% Ut for 25 cycle<br><5% Ut for 5 sec. | The quality of supply network voltage should be typical of commercial or hospital environments. If the user of the lamp MP 3020 requires continuous use even without a supply network, use an uninterruptible power supply.  |
| Power frequency magnetic field<br>IEC/EN61000-4-8  | 3A/m   | Level of magnetic field at the network frequency typical of commercial or hospital environments.   |
| Conducted immunity<br>IEC/EN61000-4-6  | 3Vrms 150kHz to 80MHz<br>(for non life-supporting equipment)                             | Portable and mobile RF communication devices should not be used near any part of the dental unit, including cables, unless they comply with recommended distances calculated with the applicable equation for transmitter frequency.<br>Recommended distances:<br>d = 1.2√P<br>d = 1.2√P from 80 Mhz to 800 MHz<br>d = 2.3√P from 800 Mhz to 2.5 GHz<br>P is the maximum nominal power issued by the transmitter in Watts (W) depending on the manufacturer of the transmitter, and d is the recommended distance in metres (m).<br>The intensity of the fixed RF transmitter field, as established in an electromagnetic investigation of site a, could be less than the compliance level of each frequency interval.<br>There can be interference near devices marked with the following symbol:  |
| Conducted immunity<br>IEC/EN61000-4-6  | 3Vrms 80MHz to 2.5GHz<br>(for non life-supporting equipment)                             |  |

Note: Ut is the power-line voltage

Note 1: The highest frequency interval is applied at 80 MHz and 800MHz.

Note 2: These guidelines might not apply to all situations. Electromagnetic propagation is influenced by absorption and reflection of structures, objects and persons.

a) ISN bands (industrial, scientific and medical) between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz and 40.66 MHz to 40.70 MHz.

b) Compliance levels in ISN bands between 150 kHz and 80 MHz and 80 MHz to 2.5 GHz present a decreasing probability of portable transmission devices causing interference if inadvertently taken to the patient area.

Therefore, an additional 10/3 factor has been incorporated into the formula used to calculate the distance between transmitters.

c) Field intensities for fixed transmitters such as base stations for radiotelephones (mobiles and cordless) and cellular mobile radios on land, CB user equipment, AM and FM transmitters and TV transmitters cannot be theoretically estimated with precision. To establish an electromagnetic environment caused by fixed RF transmitters, an electromagnetic investigation of the site should be considered. If field intensity measured at the site of use of the dental unit exceeds the aforementioned applicable compliance level, normal function of the lamp should be monitored. If any abnormal performance is noticed, additional provisions such as a different orientation or position of the lamp might be necessary.

d) The field intensity in an interval of frequencies from 150 kHz to 80 MHz should be less than 3 V/m.

**Recommended distances between portable and mobile radiocommunication devices and the dental unit**

The lamp MP 3020 is designed to function in an electromagnetic environment in which irradiating RF disturbances are under control. The client or operator of the unit can contribute toward preventing electromagnetic interferences by ensuring a minimum distance between mobile and portable RF communication devices (transmitters) and the dental unit, as recommended below, depending on the maximum output power of the radiocommunication devices.

| Maximum nominal output power of the transmitter<br>W | Distance for transmitter frequencies (m) |   |  |
|--|--|---|--|
|  | 150 kHz to 80 MHz<br>$d = 1,2 \sqrt{P}$  | 80 MHz to 800 MHz<br>$d = 1,2 \sqrt{P}$ | 800 MHz to 2,5 GHz<br>$d = 2,3 \sqrt{P}$ |
| 0,01   | 0,12                                     | 0,12                                    | 0,23                                     |
| 0,1  | 0,38                                     | 0,38                                    | 0,73                                     |
| 1  | 1,2                                      | 1,2                                     | 2,3                                      |
| 10   | 3,8                                      | 3,8                                     | 7,3                                      |
| 100  | 12                                       | 12                                      | 23                                       |

For transmitters whose maximum nominal power is not listed above, the recommended distance  $d$  in metres (m) can be calculated by using the applicable equation for the transmitter frequency, with  $P$  as maximum nominal output of the transmitter in Watts (W), depending on the manufacturer.

Notes:

The highest frequency interval is applied at 80 MHz and 800 Mhz.

These guidelines might not apply to all situations. Electromagnetic propagation is influenced by absorption and reflection of structures, objects and persons.

## 2. DESCRIPTION OF THE DEVICE

### 2.1 Description of functioning

The Light Source is made up by two LEDs.



The light is reflected by two mirrors that allow to obtain a spot of regular and uniform luminance at all frequencies of the visible spectrum and to distribute uniformly the light in the operating field without the creation of shadows or dark areas.

The adjustment of the light intensity can be done with a Joystick. The adjustment is continuous in all the range of luminance.

Simplified Maintenance thanks to the application of new technologies that keep track of the various needs in the areas of safety, ergonomics and hygiene.

The handles of the lamp are detachable to allow their sterilisation.

### 2.2 Description of the parts

- 1 - Switch
- 2 - LED
- 3 - Mirror
- 4 - Handle
- 5 - Fuses
- 6 - Transformer

Fig. 1

### 3. INSTALLATION AND CONNECTIONS

Make sure the packaging contains the following components:

- Lamp (in the version requested)
- Bag containing the lever of the switch and the wrench
- Operating Manual

The device must be cleaned before use.

The device must be installed only by qualified personnel, competent in installation of electrical equipment in medical rooms.

The power supply must always be disconnected during installation.

#### 3.1 Assembly of the lamp on the dental unit

Install the lamp by inserting the terminal pin into the specific housing of the dental unit.

Dimension of the terminal pin: 35 mm

Maximum load allowed on the head: 1 kg.

#### 3.2 Power Connections

**POWER CABLE:** In the version without transformer be sure to keep no more than the bare minimum length for the connection.

Connect the cables of the lamp to the supply with low voltage alternate current. Take care to respect the electrical features stamped on the name plate of the lamp.

#### 3.3 Remote cable connection (optional)

The connection of the remote cables must be done according with fig. 2.

Maximum range of arm from the side of the pin: 2.5 m.

The remote cable must not be lengthened during the installation phase.


Fig. 2

The connection of external devices to the remote cable could affect the compliance of the lamp to EMC requirements.

## 4. INSTRUCTIONS FOR USE

### 4.1 Switching on/off adjustment

To turn the lamp on and off, press and release the joystick lever to the left or right side.

### 4.2 Adjustment

a) to reduce light intensity, keep the joystick lever pressed to the left side (rear view of lamp) until it reaches the desired intensity. When it reaches the minimum intensity, the machine will beep once.

b) to raise the light intensity, keep the joystick lever pressed to the right side (rear view of the lamp) until it reaches the desired intensity. When it reaches the maximum intensity, the machine will beep twice.

c) to go to the minimum intensity immediately, press and release the joystick lever to the front or rear. The next time the joystick is pressed and released to the front or rear, it will revert back to the previously memorized intensity.

The lamp memorises the last setting. Each time the lamp is switched on the light intensity will be the same as when it was turned off.

The control lever must be managed with care to prevent breakage.

Do not use the control lever for the movement of the whole lamp.

### 4.3 Lamp With Remote Control (optional)

To switch the lamp on or off, press and release the button marked "A" (fig 2).

Adjustment:

a) to reduce light intensity keep pushed the button "A" until the desired intensity is reached. When the minimum intensity is reached, an acoustic signal will be heard (1 beep).

b) To increase the light intensity keep pushed the button "A" until the desired intensity is reached. When the maximum intensity is reached, an acoustic signal will be heard (2 beeps).

c) to reach the minimum light intensity immediately, press the button marked "B" (fig 2). When the minimum intensity is reached an acoustic signal will be heard (1 beep). Press the button again to restore previously selected light intensity.

The lamp memorises the last adjustment. Each time the lamp is turned on, the light intensity will be the same as when it was turned off.

#### 4.4 Overall movements

All movement of the lamp must be done through the handles of the head.  
The arm of the lamp can be moved according with following schemes:


Fig 3

Do not move the lamp or the head with the switch lever.  
Do not move the lamp or the head outside the limits described above.

#### 4.5 Acoustic Warning Signals

OTP\* = Beep 30 seconds

MAX = 1 beep - 2 beep with proximity

MIN = 1 beep

1 BEEP = On the commands

1 BEEP = At start-up

OpL\*\* = Beep 30 seconds

\* OTP: Over-temperature protection

\*\* OpL: LED charger unplugged

### 5. MAINTENANCE AND PERIODIC CHECKS

The maintenance of the lamp is made up by following operations:

- Check the absence of space between joint components of the arms (yearly)
- Check readability of data on the information plate (yearly)
- Electrical safety check according to EN 62353 (every two years)
  1. Insulation resistance
  2. leakage current
- Light checks (every five years or 10,000 hours of function)
  1. Maximum lighting intensity >35000 lux
  2. CRI decay: <20% for data on the information plate.
  3. Expected value of Blue light on spectrum emitted, measured in w/m2: <100.

## 6. CLEANING, DISINFECTION AND STERILIZATION

### 6.1 Chemical products not allowed

For the cleaning DO NOT USE products containing:

- Ammonium Hydroxide;
- Sodium Hydroxide;
- Methylene chloride;
- Methilalcohol;
- Abrasive substances;
- Trichloroethylene;
- Benzene;
- Turpentine.

Failure to comply with this indication might cause:

- risk of breakage and damage of plastic parts;
- damage of protective paints.

In case of doubt contact Faro's Service.

### 6.2 Cleaning the mirrors

It is prescribed the use these products:

- Dürr FD 366 sensitive
- FARO's Perflex

Different products can damage the mirrors

In case of doubt contact Customer Service.

Slight halos do not compromise the quality of light.

Do not use detergents containing surfactants or water repellent agents because build-up can leave halos on the mirrors and on the screen.

Do not spray any product directly on the rear of head of the lamp.

### 6.3 Cleaning of handles, arms and joints

Follow the instruction of the manufacturer of following products:

- Dürr FD 300
- Dürr FD 312
- Dürr FD 322
- Dürr FD 333
- Dürr FD 366 sensitive
- Schülke Mikrozyd Sensitive
- Saniswiss Biosanitizer
- OmniDent OmniBiozyd

Use as soft cloth for the cleaning of the parts. Remove the detergent always with a soft cloth. It is strictly forbidden to use abrasive substances, detergents with trichloroethylene, benzene, oil of turpentine or similar on any parts of the light. It is strictly forbidden to clean handles or other parts by ultrasonic bath.

### 6.4 Sterilization Of The Handles

To remove the handle, unscrew button and slip it off. To insert it, push firmly and screw on button. Handles are not provided sterile and must, therefore, be sterilised before use. Sterilize the device with standard cycles (SHPGRP) at 121°/134° C using a class B sterilizing unit. A maximum of 200 sterilization cycles is allowed.

## 7. TROUBLESHOOTING

- The lamp does not switch on** - Check on the connected of power supply.  
Check on the status of the fuses.  
If none of these causes is the case contact Faro's Service.
- Light intensity is considerably low** - Clean the screen.  
If the light intensity does not return to its original intensity contact Faro's Service.
- Handles will not come off or are difficult to remove** - Check the position of the locking screws on the handle. They should be completely open.

## 8. TECHNICAL SPECIFICATION

|                                       |                                      |
|---------------------------------------|--------------------------------------|
| Power supply voltage                  | : 17-24V ~ (22-35 V) ±10% - 50/60 Hz |
| Power absorption                      | : 23 VA                              |
| Protection against electrical hazards | : Class II                           |
| Weight (without transformer)          | : 5 kg                               |

### 8.1 Optical performances at 1000 mm from the target

|                                 |                           |
|---------------------------------|---------------------------|
| Dimensions of the spot of light | : 300 mm x 200 mm (*)     |
| Lux                             | : 3.000 to 35.000 lux (*) |
| Color Temperature               | : 5.000 K (*)             |

(\*) Typical value

The distance between the lamp and the target affect the optical performance.  
The lamp is designed to have the best performance at a distance of 1000 mm.

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## **GUARANTEE CERTIFICATE**

FARO offers the final customer a **24 month** guarantee starting from the date of purchase.

Repairs under guarantee must be performed at FARO; expenses and transport risks are at the risk of the purchaser. **Repair under guarantee is considered valid only when:**

- **all sections of the certificate have been filled in and sent in advance to FARO by Fax (039.6010540).**

The guarantee covers faults due to the bad quality of the material or manufacturing defects; in the case of valid claims, the guarantee covers free repair or replacement. **Claims for damages and/or interest are excluded.** The guarantee is not considered valid, at the sole discretion of FARO, if the fault is due to tampering, damage, incorrect use, improper maintenance and normal wear and tear.

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**FARO S.p.A.**
FARO SpA si riserva il diritto di modificare, senza preavviso, le caratteristiche indicate nel presente manuale.
FARO SpA reserves the right to change the specifications of this equipment without notice.
FARO SpA se reserve le droit de modifier, sans préavis, les caractéristiques dans ce manuel.
FARO SpA behält sich recht vor, jederzeit stillschweigend technische oder bauliche Änderung vorzunehmen.
FARO SpA se reserva el derecho de modificar sin aviso previo la características incluidas en el presente manual de uso.